

SCRIPPS NEWS

VOLUME 8 NUMBER 1

SUMMER 1994

A publication of SCRIPPS LABORATORIES, INC.

IN THIS ISSUE:

PSA-ACT Complex

This new marker for prostate cancer exhibits what serum assays for PSA lack: the ability to differentiate prostate cancer from benign prostatic hyperplasia.

Assays for PSA- α_1 -Antichymotrypsin Complex May Soon Stand Alone in the Diagnosis of Prostate Cancer

PSA-ACT

PSA assays fail to differentiate prostate cancer from benign prostatic hyperplasia

Since 1981, the number of reported prostate cancer cases in the US has increased by 50%, and the number of related deaths by 40%. Approximately 130,000 cases are reported each year, resulting in over 30,000 annual deaths.¹ As the population ages and average life-span lengthens, these numbers are expected to increase. As such, early diagnosis and effective treatment of prostate cancer are increasingly important.

To date, Prostate Specific Antigen (PSA) is the most effective marker for prostate cancer. Serum assays for PSA are used for early detection of prostate cancer, monitoring disease progress, evaluation of therapy, and detecting early recurrence of disease after radical prostatectomy.^{2,3}

The goal of PSA serum assays is to detect prostate cancer prior to metastasis, that is, while still in the organ-confined state. Detecting prostate cancer localized to the prostate gland is paramount for reliable surgical excision. PSA's ability to do this has recently come into question. An excellent review of serum assays for PSA is presented by JT Wu, PhD³; in it, PSA's short-comings are discussed and solutions proposed.

What has distinguished PSA from other cancer markers is its tissue-specificity; PSA is produced only in the prostate epithelium, making it arguably the best serum marker available today for any type of cancer.³ Although PSA is tissue-specific, it is not, however, cancer-specific. Serum levels of PSA are elevated in both prostate cancer and benign prostatic hyperplasia (BPH).^{4,5} And

frequently, only extreme elevations of PSA (>50 μ g/L) are specific for prostate cancer.⁶ As a result, difficulty arises when distinguishing prostate cancer from BPH in patients with serum PSA levels in the 4-20 μ g/L range. PSA's lack of specificity for prostate cancer can cause the BPH patient unnecessary anxiety and even unnecessary surgery.

One of the solutions to PSA's low specificity for prostate cancer, as presented by Dr. Wu³, is a relatively new discovery. Recent studies suggest the existence of a potentially better serum marker for prostate cancer, PSA- α_1 -Antichymotrypsin Complex (PSA-ACT). In a report published in November 1993, measuring serum levels of PSA-ACT improved the ability of PSA measurements to differentiate prostate cancer from BPH. Assay specificity increased from 55% to 73%, at a sensitivity of 90%. More importantly, specificity in the critical 4-20 μ g/L range was also significantly increased, with no loss of sensitivity.⁷

PSA is a 33-34 kD glycoprotein with chymotrypsin-like protease activity.⁸ Normally, very low concentrations of PSA are released into the bloodstream, approximately 0-4ng/mL. This enzymatically active form of PSA complexes with the serum protease inhibitor α_1 -antichymotrypsin (ACT) and represents the predominant form of PSA in serum.^{6,9} PSA also forms a complex with the enzyme inhibitor α_2 -macroglobulin (A₂M) in serum, but to a much lesser degree.

In addition, a significant percentage of PSA is present unbound to any inhibitor in serum. This form of PSA possesses no enzymatic activity and, thus, does not interact with serum protease inhibitors like ACT and A₂M. Recently, the different complexes of PSA with various serum protease inhibitors have been identified and characterized.^{9,10,11} Although the reported percentage of each PSA-complex varies from study to study, a table presented by EP Diamandis, MD PhD¹² nicely summarizes the various forms of PSA (see Table 1).

These studies led to the discovery that the ratio of PSA-ACT to total PSA increases in the serum of prostate cancer patients and that the percentage of free, unbound PSA decreases in prostate cancer. That these ratios remain normal in the serum of BPH patients is of particular clinical importance. In a study by Stenman, *et al*¹⁰, prostate cancer patients had a significantly higher proportion of PSA-ACT than those patients with BPH. In a related experiment, Christensson, *et al*⁶ measured the proportion of free PSA to total PSA and found the ratio to be 36% lower in prostate cancer than in BPH. Although the mechanism is not yet understood, both studies indicate that prostate malignancies release more enzymatically active PSA than other conditions. This leads to a higher percentage of PSA-ACT in the serum of prostate cancer patients.

At present, commercially available PSA assays do not consistently recognize PSA-ACT, nor do they give compatible PSA results from the same patient sample. Two excellent reviews summarize these issues, which will not be discussed in this review.^{3,13}

Two interesting reports were published in 1993 describing newly developed assays for PSA-ACT. One of the studies details a monoclonal-polyclonal antibody immunofluorometric assay that equally recognizes both free PSA and PSA-ACT.¹⁴ This mono:poly assay

has a very low detection limit of 0.002µg/L; this should prove useful in the detection of residual disease after radical prostatectomy.

The second report also describes a mono:poly immunofluorometric assay.¹⁵ What makes this assay unusual, however, is that it concurrently measures free PSA and PSA-ACT, providing separate serum levels for each analyte. This greatly simplifies the determination of the percentage of PSA-ACT to total PSA in the serum of prostate cancer patients. In fact, the authors found this assay to reduce the frequency of false-positive results due to BPH by 38%.

An issue that remains to be resolved is the calibration of the PSA-ACT assay. An assay standard containing the appropriate ratio of PSA-ACT to free PSA is inherently complex. One must take into account the normal serum ratio of PSA-ACT to free PSA, as well as any interference presented by other PSA-inhibitors like A₂M. Zhou, *et al*¹³ describe yet another sensitive mono/poly immunoassay for PSA-ACT that employs calibrators prepared by adding purified PSA to female serum that had been treated with immobilized antibody to A₂M. Using such a calibrator, they reported a recovery of >98% of PSA values throughout the standard curve. In addition, the authors cite a December 1992 international standardization conference where it was recommended that PSA standards should include serum protease inhibitors so that patient sera is more accurately reproduced.

Given the evidence presented above, it seems undeniable that the determination of serum levels of PSA-ACT is invaluable in the differential diagnosis of prostate cancer from BPH. Most importantly, PSA-ACT determinations have significantly greater specificity for prostate cancer than serum PSA evaluations alone. Optimistically, further study will show that PSA-ACT immunoassays, coupled

with other diagnostic procedures, like digital rectal exam and transrectal ultrasonography, will considerably reduce the morbidity and mortality of prostate cancer.

References

1. Crawford ED, DeAntoni EP. PSA as a Screening Test for Prostate Cancer. *Urol Clin of N Amer*, 1993, 20(4): 637-646.
2. Bjartel A, Björk T, Matikainen MT, *et al*. Production of α_1 -Antichymotrypsin by PSA-Containing Cells of Human Prostate Epithelium. *Urology*, 1993, 42(5): 502-510.
3. Wu JT. Assay for Prostate Specific Antigen (PSA): Problems and Possible Solutions. *J Clin Lab Anal*, 1994, 8: 51-62.
4. Drago JR, Badalament RA, Wientjes MG, *et al*. Relative Value of Prostate-Specific Antigen and Prostatic Acid Phosphatase in Diagnosis and Management of Adenocarcinoma of Prostate. *Urology*, 1989, 34: 187.
5. Catalona WJ, Smith DS, Ratlief TL, *et al*. Measurement of Prostate-Specific Antigen in Serum as a Screening Test for Prostate Cancer. *NE J Med*, 1991, 324: 1156.
6. Christensson A, Björk T, Nilsson O, *et al*. Serum Prostate Specific Antigen Complexed to α_1 -Antichymotrypsin as an Indicator of Prostate Cancer. *J Urol*, 1993, 150: 100-105.
7. Lilja H. Significance of Different Molecular Forms of Serum PSA. *Urol Clin of N Amer*, 1993, 20(4): 681-686.
8. Loo R, Watt KWK. *Clinical Aspects of Prostate Cancer*. Elsevier Science Publishing Co, New York, 1989: 147.
9. Lilja H, Christensson A, Dahlen U, *et al*. Prostate-Specific Antigen in Serum Occurs Predominantly in Complex with α_1 -Antichymotrypsin. *Clin Chem*, 1991, 37(9): 1618-1625.
10. Stenman UH, Leinonen J, Alfthan H, *et al*. A Complex Between Prostate-Specific Antigen and α_1 -Antichymo-trypsin is the Major Form of Prostate-Specific Antigen in Serum of Patients with Prostatic Cancer: Assay of the Complex Improves Clinical Sensitivity for Cancer. *Can Res*, 1991, 51: 222-226.
11. Wood WG, van der Sloot E, Böhle A. The Establishment and Evaluation of Luminescent-Labelled Immunometric Assays for Prostate-Specific Antigen- α_1 -Antichymotrypsin Complexes in Serum. *Eur J Clin Chem Clin Biochem*, 1991, 29: 787.
12. Diamandis EP. An Update on Prostate Specific Antigen. *Clin Chem News*, 1993, Oct: 18-19.
13. Zhou AM, Tewari PC, Bluestein BI, *et al*. Multiple Forms of Prostate-Specific Antigen in Serum: Differences in Immunorecognition by Monoclonal and Polyclonal Assays. *Clin Chem*, 1993, 39(12): 2483-2491.
14. Yu H, Diamandis EP. Ultra-sensitive Time-Resolved Immuno-fluorometric Assay of Prostate-Specific Antigen in Serum and Preliminary Clinical Studies. *Clin Chem*, 1993, 39(10): 2109-2114.
15. Leinonen J, Lövgren T, Vornanen T, *et al*. Double-Label Time-Resolved Immunofluorometric Assay of Prostate-Specific Antigen and of its Complex with α_1 -Antichymotrypsin. *Clin Chem*, 1993, 39(10): 2099-2103

Table 1. Molecular Forms of PSA in Human Serum¹²

<u>Molecular Form</u>	<u>Approximate Mol. Weight</u>	<u>% of Total</u>	<u>Comments</u>
Free PSA	30 kD	10-40	Percentage of free PSA decreases in the serum of cancer patients as total PSA increases
PSA-ACT	100 kD	60-90	Percentage of PSA-ACT complex increases as total PSA increases
PSA- α_1 -proteinase inhibitor	190 & 80 kD	< 1	Present only when PSA exceeds 40 µg/L
PSA-A ₂ M	800 kD	< 0.1	Not measured by commercial PSA kits
PSA-inter- α_1 -trypsin inhibitor	250 kD	< 0.1	Not measured by commercial PSA kits