

Metabolic Indicators

Liver Function

Catalog No.	Product	Source	Purity
B0114	Bilirubin Conjugate	Synthetic	≥97.5%

LIVER FUNCTION

Bilirubin Conjugate from Scripps is a highly-purified synthetic bilirubin derivative. It behaves identically to native bilirubin glucuronides within experimental limits. It is supplied as a lyophilized powder, blanketed with argon.

Vitamin Absorption

Catalog No.	Product	Source	Purity	Activity
F0524	Folate Binding Protein	Bovine	Reported	Reported
I1024	Intrinsic Factor	Porcine	Reported	Reported

VITAMIN ABSORPTION

This preparation of Folate Binding Protein is purified from whey protein isolated from bovine milk. It is provided as a liquid in a phosphate buffered solution. Purity, protein content, and activity are reported upon assay.

Intrinsic Factor is purified from porcine stomach and is supplied as a liquid in distilled water. Purity, protein content, and activity are reported upon assay.