



PRODUCT SPECIFICATION

Albumin, Reagent Grade, Human

Catalog Number: A0213 **Part Number:** 90076

Source: Human Serum or Plasma

Form

Lyophilized powder, essentially salt-free

Purity Specification

≥ 95%

Protein / Content Specification

Report protein to dry weight

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C frozen below -15°C. Keep material dry in a tightly sealed container. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative/non-reactive by FDA licensed tests for the following: antibodies to HIV 1 & 2, antibody to HCV, HBsAg, HIV antigen, HIV 1 RNA, and HCV RNA. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.