



PRODUCT SPECIFICATION

Avidin, Chicken Egg

Catalog Number: A1514 **Part Number:** 90081

Source: Chicken Egg White

Form

Lyophilized, essentially salt-free. Soluble in distilled water or buffer.

Purity Specification

≥ 90% protein

Activity Specification

≥ 10 Worthington Units / mg protein

One Worthington unit will bind one microgram of D-Biotin at pH 8.9.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.