



## PRODUCT SPECIFICATION

### Bilirubin Conjugate

**Catalog Number:** B0114      **Part Number:** 90604

**Source:** Synthetic Ditaurobilirubin, Disodium Salt, Molecular Formula:  $C_{37}H_{44}N_6O_{10}S_2Na_2$

#### **Form**

Lyophilized orange to reddish brown powder.

#### **Activity Specification**

Molar Absorptivity @ 450nm in Water:  $\geq 55,000$  (anhydrous salt) corrected for water and residual solvent.

(anhydrous MW = 842.9)

#### **Characteristics Specification**

Report UV-VIS Spectrum.

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store frozen below  $-20^{\circ}C$ . Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*