



## **PRODUCT SPECIFICATION**

### **Conjugate Dilution Buffer**

**Catalog Number:** B0221      **Part Number:** 90092

#### **Form**

Clear yellow-golden liquid; buffer containing proprietary stabilizing ingredients.

#### **Characteristics Specification**

This product is designed for the resuspension and dilution of colloidal gold conjugates. The buffer promotes the stabilization of colloidal gold conjugates dried on the conjugate matrix, and provides for rapid reconstitution.

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures. The buffer should not contact eyes, skin, or clothing. Wear gloves, lab coat, and eye protection when handling. Wash hands thoroughly after use.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*