



PRODUCT SPECIFICATION

Carcinoembryonic Antigen (CEA), Human, $\geq 50\%$

Catalog Number: C0223 **Part Number:** 90097

Source: Human liver metastasis of a primary colon carcinoma

Form

Clear liquid - Buffer: 10mM Tris, 0.1% Sodium Azide, pH 8.0 \pm 0.1

Purity Specification

$\geq 50\%$ (based on ratio of specific protein by immunoassay to total protein by Lowry Protein Assay)

Protein / Content Specification

Report specific protein by immunoassay (specify method). Report total protein by Lowry Protein Assay (BSA standard).

Product is packaged based on the protein concentration by immunoassay.

Characteristics Specification

Indicate calibration standard for immunoassay, if available

Storage / Shipping

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.