



PRODUCT SPECIFICATION

Ferritin, Human Spleen

Catalog Number: F0424 **Part Number:** 90007

Source: Human Spleen

Form

Amber to red liquid - Buffer: 10mM Tris, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 8.0 ± 0.1

Purity Specification

≥ 98% by SDS-PAGE

Protein / Content Specification

≥ 2 mg/ml by immunoassay (report method)

≥ 2 mg/ml by Lowry Protein Assay (BSA standard)

Product is packaged based on the immunoassay value

Characteristics Specification

Indicate calibration standard for the immunoassay (if available)

Storage / Shipping

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.