



## PRODUCT SPECIFICATION

### Follicle Stimulating Hormone (FSH), Human, $\beta$ -Subunit

**Catalog Number:** F0814      **Part Number:** 90107

**Source:** Human Pituitary Glands

**Form**

White to off-white powder, lyophilized from 50mM ammonium bicarbonate

**Purity Specification**

≥ 98% by SDS-PAGE

**Protein / Content Specification**

Determined prior to lyophilization

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*