



PRODUCT SPECIFICATION

Goat antibody to Human Creatine Kinase MM (CKMM), purified

Catalog Number: GC139 **Part Number:** 90123

Source: Goat Plasma

Immunogen: Purified CKMM from Human Skeletal Muscle

Purification Method: Affinity chromatography on immobilized human CKMM

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 500mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 ± 0.1

Purity Specification

Report purity by SDS-PAGE

Protein / Content Specification

≥ 5 mg/ml by A₂₈₀ (E_{mg/ml} = 1.38)

Characteristics Specification

Report IEP results vs. CKMM, CKMB, CKBB and 2x human plasma

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.