

PRODUCT SPECIFICATION

Goat antibody to Human Ferritin, IgG fraction

Catalog Number: GF035 Part Number: 90560

Source: Goat Serum

Immunogen: Purified Ferritin from Human Liver

Purification Method: Precipitation

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.09% Sodium Azide,

 $pH 7.2 \pm 0.1$

Protein / Content Specification

Report total protein concentration by A280

Activity Specification

Report IEP results vs. Human Liver Report IEP results vs. 2x Human Plasma

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.

Tel: 858.546.5800 Fax: 858.546.5812