



## PRODUCT SPECIFICATION

### Horseradish Peroxidase (HRP) Conjugate Stabilizer, no preservative

**Catalog Number:** H1521      **Part Number:** 90135

#### **Form**

Yellow to tan liquid. Preparation is a phosphate buffered saline solution containing bovine serum albumin and other proprietary ingredients. Product has been 0.2 $\mu$ m filtered into gamma-irradiated bottles. Buffer components may crystallize at low temperature. Warm in a 30°C water bath and invert several times until crystals dissolve.

#### **Characteristics Specification**

HRP-Conjugate Stabilizer is a ready-to-use product designed for the reconstitution of lyophilized HRP conjugates. Conjugates prepared in this manner may retain peroxidase activity for up to 2 years when stored at 2°C to 8°C. This product does not contain a preservative. It must be handled using aseptic techniques to ensure proper performance.

Expiration Date:      (report 18 months from manufacturing date)

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*