



PRODUCT SPECIFICATION

Horseradish Peroxidase (HRP) Conjugate Stabilizer w/proprietary preservative

Catalog Number: H1522 **Part Number:** 90136

Form

Yellow to tan liquid. Preparation is a phosphate buffered saline solution containing bovine serum albumin and other proprietary ingredients, including a proprietary preservative. Product has been 0.2µm filtered into gamma-irradiated bottles. Buffer components may crystallize at low temperature. Warm in a 30°C water bath and invert several times until crystals dissolve.

Characteristics Specification

HRP-Conjugate Stabilizer is a ready-to-use product designed for the reconstitution of lyophilized HRP conjugates. Conjugates prepared in this manner may retain peroxidase activity for up to 2 years when stored at 2°C to 8°C. Although a preservative has been added to this product, it should be handled using aseptic techniques to ensure proper performance.

Expiration Date: (report 18 months from manufacturing date)

Storage / Shipping

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.