



PRODUCT SPECIFICATION

Human Immunoglobulin E (IgE) Myeloma Serum, Kappa

Catalog Number: I0223 **Part Number:** 90139

Source: Human Serum (Myeloma)

Form

Liquid, neat serum without preservative

Protein / Content Specification

Report IgE concentration

Activity Specification

Report activity (IU/ml) and method

Characteristics Specification

Indicate calibration standard for assay, if available

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.