



PRODUCT SPECIFICATION

Human Immunoglobulin E (IgE lambda), Plasma

Catalog Number: I0324 **Part Number:** 90597

Source: Normal Human Plasma spiked with Cell Culture Human IgE, lambda

Form

Liquid, neat human plasma (EDTA) prepared by addition of human IgE lambda (cell culture). Product does not contain preservative.

Protein / Content Specification

>15,000 IU/ml (report method)
Report IgE concentration (2.4 ng/IU)

Characteristics Specification

Indicate calibration standard for assay, if available

Storage / Shipping

Product is shipped using dry ice (***not temperature monitored***). On arrival, store frozen below -10°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.