



## PRODUCT SPECIFICATION

### Immunoglobulin G (IgG), Goat

**Catalog Number:** I0514      **Part Number:** 90142

**Source:** Goat Serum

**Form**

White to off-white lyophilized powder

**Purity Specification**

Report purity

**Protein / Content Specification**

Report approximate % protein by  $A_{280}$  / dry weight  
Balance is buffer salts.

**Characteristics Specification**

This product is bulk processed and packaged by dry weight. Protein is determined by  $A_{280}$  ( $E_{mg/ml} = 1.3$ ).

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*