



PRODUCT SPECIFICATION

Immunoglobulin G (IgG), Rabbit, Lyophilized - Gamma Globulin Fraction

Catalog Number: I0713 **Part Number:** 90587

Source: Rabbit Serum

Form

White to off-white powder

Purity Specification

≥ 95% by Serum Protein Electrophoresis

Protein / Content Specification

≥ 85% protein by A_{280} / dry weight

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.