



PRODUCT SPECIFICATION

Intrinsic Factor, porcine

Catalog Number: I1024 **Part Number:** 90040

Source: Porcine Stomach

Form

Clear liquid, supplied in distilled water. Dilute with 100mM Potassium Phosphate, pH 7.5 ± 0.1, containing 0.1mg/ml albumin

Purity Specification

Report purity and method

Protein / Content Specification

Report protein concentration

Activity Specification

Report B12 binding ability/mg protein

One unit will bind 1ng of Vitamin B12

Storage / Shipping

Product is shipped using dry ice (*not temperature monitored*). On arrival, store frozen below -20°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.