



PRODUCT SPECIFICATION

Epidermal Keratin, Human

Catalog Number: K0124 **Part Number:** 90151

Source: Human Epidermis

Form

Clear liquid - Buffer: 50mM Tris, 8M Urea, 100mM β -Mercaptoethanol, 0.1% Sodium Azide, pH 8.4 \pm 0.1

Purity Specification

\geq 90%

Protein / Content Specification

Report protein concentration
(modified Lowry Protein Assay, BSA standard)

Storage / Shipping

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

Tested and found negative for HIV 1 and Hepatitis B. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.