



PRODUCT SPECIFICATION

β -2-Microglobulin (β 2M), Human, \geq 40%

Catalog Number: M0113 **Part Number:** 90160

Source: Human urine from patients with chronic renal tubular proteinuria

Form

White to off-white powder, lyophilized from 20mM ammonium bicarbonate

Purity Specification

\geq 40% (specify method)

Protein / Content Specification

Determined by RID prior to lyophilization

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.