



## PRODUCT SPECIFICATION

### Monoclonal antibody to Human C-Reactive Protein (CRP)

**Catalog Number:** MC017      **Part Number:** 90274

**Clone Number:** CC002

**Source:** Murine

**Immunogen:** Purified human C-Reactive Protein

**Purification Method:** Protein A Chromatography and gel filtration

#### **Form**

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.4 ± 0.1

#### **Purity Specification**

Protein A purified IgG fraction.  
Report purity by SDS-PAGE

#### **Protein / Content Specification**

≥ 5 mg/ml by A<sub>280</sub> (E<sub>mg/ml</sub> = 1.4)

#### **Characteristics Specification**

Mouse IgG<sub>2b,λ</sub> isotype

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*