



PRODUCT SPECIFICATION

Monoclonal antibody to Human Creatine Kinase BB (CKBB)

Catalog Number: MC117 **Part Number:** 90177

Clone Number: BC006

Source: Murine

Immunogen: Purified human CKBB

Purification Method: Protein A Chromatography

Form

Clear liquid - Buffer: Phosphate Buffered Saline, 0.05% Sodium Azide, pH 7.2 ± 0.1

Purity Specification

Protein A purified IgG fraction.
Report purity by SDS-PAGE

Protein / Content Specification

Report protein concentration by A_{280} ($E_{mg/ml} = 1.4$)

Characteristics Specification

Mouse IgG_{2a,k}

This clone recognizes the CK-B subunit in both CK-BB and CK-MB.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.