



PRODUCT SPECIFICATION

Monoclonal antibody to Human Ferritin

Catalog Number: MF037 **Part Number:** 90582

Clone Number: 201601

Source: Murine

Immunogen: Purified human ferritin

Purification Method: Protein A Chromatography

Form

Clear liquid - Buffer: 15mM Phosphate Buffered Saline, 0.1% Sodium Azide, pH 7.4 ± 0.1

Purity Specification

Protein A purified IgG fraction

Protein / Content Specification

Report protein concentration by A_{280} ($E_{mg/ml} = 1.4$)

Characteristics Specification

Mouse IgG₁ isotype

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.