

Product Review

Raw Materials for Immunoassay Development and Manufacturing

Animal Immunoglobulins - Reagents to Reduce Non-Specific Binding

Purified Immunoglobulins


Catalog No.	Description	Purity
I0514	Goat IgG	≥95%
I0624	Mouse IgG	≥95%
I0724	Rabbit IgG	≥95%

IMMUNOGLOBULINS

Non-specific immunoglobulins, primarily IgG, are often used as blocking agents to minimize non-specific binding in an immunoassay and thereby reduce signal background.

Scripps offers purified IgGs from various species for use in immunoassay development. The materials are purified from animal sera and are offered with consistent high quality and performance.

Example SDS-PAGE Gel of Mouse IgG (lane 5)



Lane	Sample
1	Molecular Weight Markers
3	Side Fraction of Mouse IgG Production
5	Purified Pool of Mouse IgG for Final Lot