

# Product Review

Raw Materials for Immunoassay Development and Manufacturing

## Metabolic Indicators

### Liver Function

Catalog No.	Product	Source	Purity
B0114	Bilirubin Conjugate	Synthetic	≥97.5%

### LIVER FUNCTION

Bilirubin Conjugate from Scripps is a highly-purified synthetic bilirubin derivative. It behaves identically to native bilirubin glucuronides within experimental limits. It is supplied as a lyophilized powder, blanketed with argon.

### Vitamin Absorption

Catalog No.	Product	Source	Purity	Activity
F0524	Folate Binding Protein	Bovine	Reported	Reported
I1024	Intrinsic Factor	Porcine	Reported	Reported

### VITAMIN ABSORPTION

This preparation of Folate Binding Protein is purified from whey protein isolated from bovine milk. It is provided as a liquid in a phosphate buffered solution. Purity, protein content, and activity are reported upon assay.

Intrinsic Factor is purified from porcine stomach and is supplied as a liquid in distilled water. Purity, protein content, and activity are reported upon assay.