



PRODUCT SPECIFICATION

Rabbit antibody to Human Cathepsin D, IgG fraction

Catalog Number: RC245 **Part Number:** 90236

Source: Rabbit Serum

Immunogen: Purified Cathepsin D from Human Liver

Purification Method: Ion-Exchange Chromatography

Form

Clear liquid - Buffer: 10mM Phosphate, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 ± 0.1

Purity Specification

DEAE purified IgG fraction

Protein / Content Specification

Report total protein concentration

Characteristics Specification

Report lot specific data (i.e. ELISA)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.