



## PRODUCT SPECIFICATION

### Sex Hormone Binding Globulin (SHBG), High Purity

**Catalog Number:** S1725      **Part Number:** 90242

**Source:** Human Serum

**Form**

Clear liquid - Buffer: 20mM Tris-HCl, 10% Glycerol, pH 7.4 ± 0.1

**Purity Specification**

≥ 95% by SDS-PAGE

Report contaminants

**Protein / Content Specification**

Report protein concentration by RIA and/or Lowry

**Storage / Shipping**

Product is shipped using dry ice (*not temperature monitored*). On arrival, store frozen below -60°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*