



PRODUCT SPECIFICATION

Thyroid Stimulating Hormone, β -Subunit (β -TSH), Human

Catalog Number: T0314 **Part Number:** 90247

Source: Human Pituitary Glands

Form

White to off-white powder, lyophilized from 50mM ammonium bicarbonate

Purity Specification

≥ 90% by SDS-PAGE

Report contaminant levels if available

Protein / Content Specification

Determined prior to lyophilization
(Lowry Protein Assay, BSA standard)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Source materials have not been tested for infectious diseases. This product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.