



PRODUCT SPECIFICATION

Thyroglobulin, Human, low IgG

Catalog Number: T0714 **Part Number:** 90349

Source: Human Thyroid Glands

Form

White to tan powder, lyophilized from 20mM ammonium bicarbonate.

Purity Specification

≥ 96% (report method)

Material has been adsorbed to decrease IgG levels

Protein / Content Specification

Report protein concentration prior to lyophilization by A_{280} ($E_{mg/ml} = 1.04$)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.