



## PRODUCT SPECIFICATION

### Troponin C, Human Heart

**Catalog Number:** T5224      **Part Number:** 90338

**Source:** Human Heart

**Form**

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.02% Sodium Azide, pH 7.4 ± 0.1

**Purity Specification**

≥ 95% by SDS-PAGE

**Protein / Content Specification**

≥ 0.5 mg/ml by the Bradford method  
(BioRad Protein Assay, BSA standard)

**Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*