

PRODUCT SPECIFICATION

α-Fetoprotein (AFP), Affinity Purified, Human, ≥ 95%

Catalog Number: A0724 **Part Number:** 90492

Source: Human Cord Serum

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, 0.05% Sodium Azide, pH 7.2 ± 0.1.

Purity Specification

≥ 95% by SDS-PAGE

Protein / Content Specification

≥ 1.0 mg/ml by Immunoassay.

Characteristics Specification

Indicate calibration standard for immunassay, if available

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.