



LABORATORIES

Manufacturers of the Highest Quality

PRODUCT SPECIFICATION

NT-ProBNP (N-Terminal Brain, or B-Type, Natriuretic Peptide), Recombinant

Catalog Number: B0414 **Part Number:** 90625

Source: *E.coli*

Form

White powder, lyophilized from 200mM ammonium acetate

Purity Specification

≥ 95% by SDS-PAGE

Protein / Content Specification

Report protein concentration prior to lyophilization by A₂₈₀

Storage / Shipping

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812
scrippslabs.com | contact@scrippslabs.com