



## PRODUCT SPECIFICATION

### CEA and PSA Dilution Buffer Tris/Azide, pH 8.0

**Catalog Number:** B1012      **Part Number:** 90548

#### **Form**

Supplied as a clear liquid, 0.2µm filtered.

Formulation: 10mM Tris, 0.1% Sodium Azide, pH 8.0 ± 0.1

#### **Characteristics Specification**

CEA & PSA Dilution Buffer is a ready-to-use product intended for the dilution of purified CEA (cat. # C0223 and C0224), PSA (cat. # P0723, P0724, P0725, and P0715), and the non-complexing form of PSA (cat. # P5125). After diluting, the protein should be stored aliquots as indicated on the certificate.

Although a preservative has been added to this product and the material has been filtered (0.2µm) prior to bottling in nonpyrogenic and noncytotoxic media bottles, it is not packaged and tested for sterility. The dilution buffer should be re-filtered for use in sterile applications.

#### **Storage / Shipping**

Product is shipped at ambient temperature or using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*