



## PRODUCT SPECIFICATION

### C-Reactive Protein (CRP), Human Fluids, $\geq 99\%$

**Catalog Number:** C0129      **Part Number:** 90518

**Source:** Human Fluids

**Form**

Clear liquid - Buffer: 20mM Tris, 280mM Sodium Chloride, 5mM Calcium Chloride, 0.1% Sodium Azide, pH  $8.0 \pm 0.1$

**Purity Specification**

$\geq 99\%$  by SDS-PAGE

**Protein / Content Specification**

$\geq 2.0$  mg/ml by  $A_{280}$  ( $E_{mg/ml} = 1.7$ )

**Storage / Shipping**

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at  $2^{\circ}\text{C}$  to  $8^{\circ}\text{C}$ . Do not freeze!

**Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*