



LABORATORIES

Manufacturers of the Highest Quality

PRODUCT SPECIFICATION

C-Reactive Protein (CRP), Human Fluids, $\geq 99\%$

Catalog Number: C0129 **Part Number:** 90518

Source: Human Fluids

Form

Clear liquid - Buffer: 20mM Tris, 280mM Sodium Chloride, 5mM Calcium Chloride, 0.1% Sodium Azide, pH 8.0 \pm 0.1

Purity Specification

$\geq 99\%$ by SDS-PAGE

Protein / Content Specification

≥ 2.0 mg/ml by A₂₈₀ (E_{mg/ml} = 1.7)

Storage / Shipping

Product is shipped using a refrigerated gel pack (*not temperature monitored*). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812
scrippslabs.com | contact@scrippslabs.com