

PRODUCT SPECIFICATION

Human Chorionic Gonadotropin (hCG, human), ≥ 95%, Recombinant

Catalog Number: C0733 **Part Number:** 90646

Source: Intact human hCG (α- & β-Subunits) produced in HEK 293 cells

Form

White powder, lyophilized from 50mM ammonium bicarbonate

Purity Specification

≥ 95% by SDS-PAGE

Protein / Content Specification

Report protein concentration prior to lyophilization
(Lowry Protein Assay, BSA standard)

Activity Specification

Report IU & mg protein per vial by Centaur CP & Lowry
Product is packaged by Centaur CP immunoassay value.

Target values:

650 IU / vial; 1,300 IU / vial; 13,000 IU / vial

Storage / Shipping

Product is shipped using refrigerated gel packs (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.