



## PRODUCT SPECIFICATION

### Ferritin, Human Liver, Immunoaffinity Purified

**Catalog Number:** F0325      **Part Number:** 90563

**Source:** Human Liver

#### **Form**

Amber to red liquid - Buffer: 10mM Tris, 150mM Sodium Chloride, 0.1% Sodium Azide, pH 8.0 ± 0.1

#### **Purity Specification**

≥ 98% by SDS-PAGE

#### **Protein / Content Specification**

1.0 - 5.0 mg/ml by immunoassay (specify method)

1.0 - 5.0 mg/ml by Lowry Protein Assay (BSA standard)

Product is packaged based on the immunoassay value

#### **Characteristics Specification**

Indicate calibration standard for the immunoassay (if available)

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

**Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C.** Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*