



## PRODUCT SPECIFICATION

### Folate Binding Protein

**Catalog Number:** F0524      **Part Number:** 90044

**Source:** Bovine Milk

**Form**

Clear liquid - Buffer: 120mM Sodium Phosphate, 30mM Sodium Chloride, pH 7.4 ± 0.1

**Purity Specification**

Report results and method

**Protein / Content Specification**

Report protein concentration. This is determined indirectly using the measured activity in µg/ml and the binding capacity of 8 µg folic acid/mg folate binding protein

**Activity Specification**

Report folate binding ability in µg/ml

One mg of active folate binding protein will bind 8µg of folic acid

**Storage / Shipping**

Product is shipped using dry ice (***not temperature monitored***). On arrival, store frozen below -20°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

**Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*