



## PRODUCT SPECIFICATION

### Folate Binding Protein

**Catalog Number:** F0524      **Part Number:** 90589

**Source:** Whey protein isolated from bovine milk

**Form**

Clear liquid - Buffer: 120mM Sodium Phosphate, 30mM Sodium Chloride, pH 7.4 ± 0.1

**Purity Specification**

Report purity by HPLC

**Protein / Content Specification**

Report protein concentration by Lowry Assay

**Activity Specification**

Report folate binding activity by Centaur CP (µg folic acid / mg protein)

**Storage / Shipping**

Product is shipped using dry ice (*not temperature monitored*). On arrival, store frozen below -20°C. Avoid repeated freeze-thaw cycles by packaging into smaller sizes and freezing.

**Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*