



L A B O R A T O R I E S

Manufacturers of the Highest Quality

PRODUCT SPECIFICATION

Follicle Stimulating Hormone (FSH, human), $\geq 95\%$, Recombinant

Catalog Number: F0617 **Part Number:** 90632

Source: Intact human FSH (α - & β -Subunits) produced in HEK 293 cells

Form

White powder, lyophilized from 50mM ammonium bicarbonate

Purity Specification

$\geq 95\%$ by SDS-PAGE

Protein / Content Specification

Report protein concentration prior to lyophilization
(Lowry Protein Assay, BSA standard)

Activity Specification

Report IU & mg protein per vial by Centaur CP & Lowry
Product is packaged by Centaur CP immunoassay value.

Target values:

70 IU / vial; 350 IU / vial; 700 IU / vial; 7,000 IU / vial

Storage / Shipping

Product is shipped using refrigerated gel packs (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812
scrippslabs.com | contact@scrippslabs.com