



LABORATORIES

Manufacturers of the Highest Quality

PRODUCT SPECIFICATION

Follicle Stimulating Hormone (FSH), Human, β -Subunit

Catalog Number: F0814 **Part Number:** 90107

Source: Human Pituitary Glands

Form

White to off-white powder, lyophilized from 50mM ammonium bicarbonate

Purity Specification

≥ 98% by SDS-PAGE

Protein / Content Specification

Determined prior to lyophilization

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812
scrippslabs.com | contact@scrippslabs.com