

PRODUCT SPECIFICATION

Follicle Stimulating Hormone (FSH), Porcine

Catalog Number: F0913 **Part Number:** 90108

Source: Porcine Pituitary Glands

Form

Lyophilized powder

Activity Specification

≥ 0.5 x NIH-FSH-S1

Characteristics Specification

This product is bulk processed and packaged by dry weight

Activity is based on an hCG-augmented ovarian weight-gain assay. Endocrinology, 69: 604 (1953)

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.