



## PRODUCT SPECIFICATION

### Human Apoferritin, H+L Chains, $\geq 95\%$ , Recombinant

**Catalog Number:** F1021      **Part Number:** 90659

**Source:** Human Ferritin (Heavy & Light Chains), produced in *E.coli*

#### **Form**

Clear liquid - Buffer: 10mM Tris, 150mM Sodium Chloride, 0.02% Sodium Azide, pH 8.0  $\pm$  0.1

#### **Purity Specification**

$\geq 95\%$  by SDS-PAGE

#### **Protein / Content Specification**

1.0 - 5.0 mg/ml by immunoassay (specify method)

1.0 - 5.0 mg/ml by Lowry Protein Assay (BSA standard)

Product is packaged based on the immunoassay value

#### **Characteristics Specification**

Indicate calibration standard for the immunoassay (if available)

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*