

## **PRODUCT SPECIFICATION**

### **Human Apoferritin, H+L Chains, ≥ 95%, Recombinant**

**Catalog Number:** F1021      **Part Number:** 90659

**Source:** Human Ferritin (Heavy & Light Chains), produced in *E.coli*

**Form**

Clear liquid - Buffer: 10mM Tris, 150mM Sodium Chloride, 0.02% Sodium Azide, pH 8.0 ± 0.1

**Purity Specification**

≥ 95% by SDS-PAGE

**Protein / Content Specification**

1.0 - 5.0 mg/ml by immunoassay (specify method)

1.0 - 5.0 mg/ml by Lowry Protein Assay (BSA standard)

Product is packaged based on the immunoassay value

**Characteristics Specification**

Indicate calibration standard for the immunoassay (if available)

**Storage / Shipping**

Product is shipped using a refrigerated gel pack (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C. Do not freeze!

**Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.