



PRODUCT SPECIFICATION

Growth Hormone (GH), Human, $\geq 70\%$

Catalog Number: G0714 **Part Number:** 90045

Source: Human Pituitary Glands

Form

White powder, lyophilized from 50mM ammonium bicarbonate

Purity Specification

$\geq 70\%$ by SDS-PAGE

$\leq 1\%$ FSH, LH, TSH, and Prolactin (w/w)

Protein / Content Specification

Determined by the Bradford method prior to lyophilization
(BioRad Protein Assay, BSA standard)

Activity Specification

Report activity by immunoassay (indicate method)

Characteristics Specification

Indicate calibration standard for immunoassay, if available

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

Tested and found negative for HIV 1 & 2, Hepatitis B, and Hepatitis C. Although source materials have been tested for certain diseases, no test can absolutely assure the absence of all infectious agents. Therefore, this product should be handled as a potential biohazard.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Manufacturers of the highest quality proteins and antibodies.