



## PRODUCT SPECIFICATION

### Galectin-3, Recombinant

**Catalog Number:** G1214      **Part Number:** 90624

**Source:** *E.coli* (Human Galectin-3)

#### **Form**

White powder, lyophilized from 50mM ammonium bicarbonate

#### **Purity Specification**

≥ 95% by SDS-PAGE

#### **Protein / Content Specification**

Report protein concentration prior to lyophilization by the BCA Protein Assay

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting, store frozen in small aliquots. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*