

## PRODUCT SPECIFICATION

### **Galectin3 Recombinant, 95%**

Catalog Number: G1214      Part Number: 90712

Source: *E.coli* (Human Galectin-3)

#### Form

White powder, lyophilized from 50mM ammonium bicarbonate.

Recommended resuspension buffer: 20 mM Sodium Phosphate, 150 mM Sodium Chloride, pH 7.5 ± 0.1

#### Purity Specification

**95% by SDS PAGE**

#### Protein / Content Specification

**Report protein concentration prior to lyophilization by A280**

#### Storage / Shipping

**Product is shipped using a refrigerated gel pack (not temperature monitored). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen below -15°C. After reconstituting store frozen in small aliquots. Avoid repeated freeze thaw cycles.**

#### Precaution / Hazard Statement

**This product is not considered a biohazard. Handle using normal laboratory procedures.**

Intended Use: **For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures**