



LABORATORIES

Manufacturers of the Highest Quality

PRODUCT SPECIFICATION

Goat antibody to Human Albumin, purified

Catalog Number: GA029 **Part Number:** 90116

Source: Goat Serum

Immunogen: Purified Albumin from Human Serum/Plasma

Purification Method: Affinity chromatography on immobilized human albumin

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 500mM Sodium Chloride, 0.1% Sodium Azide, pH 7.4 ± 0.1

Purity Specification

Report purity by SDS-PAGE

Protein / Content Specification

≥ 3 mg/ml by A₂₈₀ (E_{mg/ml} = 1.38)

Characteristics Specification

Report IEP results

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

6838 Flanders Drive, San Diego, CA 92121-2904 USA | +1 858.546.5800 | fax: +1 858.546.5812
scrippslabs.com | contact@scrippslabs.com