



## PRODUCT SPECIFICATION

### Goat antibody to Human C-Reactive Protein (CRP), purified

**Catalog Number:** GC019      **Part Number:** 90119

**Source:** Goat Serum

**Immunogen:** Purified CRP from Human Serum

**Purification Method:** Affinity chromatography on immobilized human CRP

#### **Form**

Clear liquid - Buffer: 10mM Sodium Phosphate, 500mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 ± 0.1

#### **Purity Specification**

Affinity purified, free from particulates.  
Report purity by SDS-PAGE.

#### **Protein / Content Specification**

≥ 5 mg/ml by A<sub>280</sub> (E<sub>mg/ml</sub> = 1.38)

#### **Characteristics Specification**

Report IEP results

#### **Storage / Shipping**

Product is shipped using a refrigerated gel pack or dry ice (***not temperature monitored***). On arrival, store refrigerated at 2°C to 8°C or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

#### **Precaution / Hazard Statement**

This product is not considered a biohazard. Handle using normal laboratory procedures.

**Intended Use:** For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

*Manufacturers of the highest quality proteins and antibodies.*