

PRODUCT SPECIFICATION

Goat antibody to Human Chorionic Gonadotropin, β -Subunit (β -hCG), purified

Catalog Number: GC099 **Part Number:** 90062

Source: Goat Serum

Immunogen: Purified Human Chorionic Gonadotropin, β -Subunit

Purification Method: Affinity chromatography on immobilized β -hCG

Form

Clear liquid - Buffer: 10mM Sodium Phosphate, 500mM Sodium Chloride, 0.1% Sodium Azide, pH 7.2 \pm 0.1

Purity Specification

Affinity purified, free from particulates.

Report purity by SDS-PAGE.

Protein / Content Specification

≥ 3 mg/ml by A280 (E_{mg/ml} = 1.38)

Characteristics Specification

Pairs with monoclonal antibody clone BC001

Report contaminants, titer, cross-reactivity, and any other batch data if available.

Storage / Shipping

Product is shipped using a refrigerated gel pack or dry ice (**not temperature monitored**). On arrival, store refrigerated at 2°C to 8°C short term (up to 2 weeks) or frozen in aliquots below -15°C. Avoid repeated freeze-thaw cycles.

Precaution / Hazard Statement

This product is not considered a biohazard. Handle using normal laboratory procedures.

Intended Use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.